Genesee

BIOMEDICAL, INC.

Design beyond standard

INSTRUCTIONS FOR USE

**Aortic Repair Tool Kit, Model ART-2,
Model ASR-19/21, Model ASR-23/25, Model ASR 27/29, Model ASR 31/33, Model ARC-6,
Model ARC-7, Model ARC-8, Model ARC-9 and Model ARC-10**

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**Instrumenten-Set für die Aortenreparatur, Modell ART-2,
Modell ASR-19/21, Modell ASR-23/25, Modell ASR 27/29, Modell ASR 31/33,
Modell ARC-6, Modell ARC-7, Modell ARC-8, Modell ARC-9 und Modell ARC-10**

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DESCRIPTION

The Aortic Repair Tool Kit consists of a set of tools with key components to evaluate the aortic valve leaflets, determine the annular size, and place and secure an annular stitch at the correct diameter of the basal ring. These specialized tools, designed for use in aortic valve remodeling and repair, include cusp calipers and annular sizer (obturators). The tools are provided in a sterilization tray for use in aortic valve or aortic root surgery.

Accurate assessment and evaluation of cusp and root geometry is an essential component for achieving functional results with aortic valve repair. Repaired aortic valves with low effective height have been shown to progress in regurgitation. In a large series of patients undergoing valve-preserving aortic root surgery, postoperative effective cusp height less than 9mm was associated with late AR and reoperation. (Bierbach et al, 2010, Kunihara et al, 2012)

The CUSP CALIPER is designed to determine intraoperative effective cusp height and geometric cusp height. The CUSP CALIPERS are provided in sizes 10, 9, 8, 7 and 6 mm for use with most patients.

The ANNULAR SIZER consists of a double-ended ANNULAR SIZER and is used to measure the aortic annulus/aorto-ventricular junction. The SIZER is used to assess the diameter of the aorto-ventricular junction. The SIZER can subsequently be used to secure a supporting stitch in the aorto-ventricular junction if applicable. The SIZER is provided in four sizes: 19/21, 23/25, 27/29 and 31/33.

PACKAGING

The CUSP CALIPERS and ANNULAR SIZERS are provided NON-STERILE and must be cleaned and sterilized prior to initial use and each subsequent use. See Cleaning and Sterilization section.

STORAGE

Store in a dry place at room temperature to protect the product and minimize contamination potential.

INDICATIONS

The CUSP CALIPERS and ANNULAR SIZERS are intended to aid in aortic valve repair or root remodeling. The tools are used to evaluate anatomical and geometric relationships in the aortic valve and aortic root, and may also be used to narrow the annulus further in aortic valve reimplantation.

METHOD OF USE

Sizing Effective Cusp Height

In every normal aortic valve there is a difference between the height of the free margins in the center of the cusp and the plane of the nadirs of the aortic insertion line. This caliper is designed to measure that difference in height; this height difference is called effective height (eH). (Figure 1, adapted from Bierbach et al.)

The caliper is placed such that the longer end rests on the lowest (i.e., central) point of the aortic insertion line. The axis of the caliper should be in line with the axis of the left ventricular outflow tract. The shorter end should then tightly rest on the free margin, with the curve accommodating

the margin. If a tight fit can be achieved with the 10 mm caliper, effective height is 10 mm; if the 8 mm caliper provides a tight fit, effective height is 8 mm. Lower values can be estimated visually in relation to the shorter end of a given size.


Figure 1

STJ: sinotubular junction; AN: annulus; eH: effective height; gH: geometric height

Procedure for Sizing Geometric Cusp Height

Normal aortic valves have been found to have a predictable amount of leaflet tissue, which is necessary to provide adequate valve function. The determination of tissue height in the central part of the cusp has been propagated as an easy indicator of amount of cusp tissue. This tissue height has been called geometric height (gH) as opposed to effective height, the height difference in the closed position. (Schäfers et al, 2012)

The commissures are placed under radial tension using stay sutures. The free margin of the cusp to be measured is held with a forceps and placed under tension. The geometric height is measured with the ruler starting from the aortic insertion in the nadir of the sinus to the nodulus of Arantius of the respective cusp. This measurement should be repeated for accuracy. (See Figure 2)


Figure 2

Assessment of Geometric Leaflet Length

Sizing the Annulus

The size of the aortic annulus has been shown to vary and enlargement has been associated with native aortic valve regurgitation as well as failure of aortic valve repair (Aicher et al, 2013, Kuniyama et al, 2013). Experience has shown that there may be marked discrepancy between

annular size measured by echocardiography and that determined by direct examination. Exact measurement is felt to be a prerequisite for standardized aortic valve repair and reproducible outcome (Lansac et al, 2010). The SIZERS are designed to allow for easy and exact intraoperative measurement of the annulus. SIZERS are entered into the aorto-ventricular junction until a snug fit is established, determining the actual annular dimension.

The SIZER can also be used in place of a Hegar dilator to determine a predicted postoperative diameter of a suture annuloplasty at the level of the basal ring. (Kunihara et al, 2012. Svensson et al, 2003) For this purpose, first the desired procedure of valve-preserving root replacement is completed. A suture is then placed at the level of the basal ring, and a SIZER of desired size is chosen. The SIZER is inserted into the aortic root and left ventricular outflow tract; the suture is then tied around the root with the SIZER in place. The SIZER is then removed. (Aicher et al, 2013)

WARNINGS

Take care not to damage the aorta, aortic valve or valve leaflets during use.

PRECAUTIONS

Take care not to damage the CUSP CALIPER or ANNULAR SIZER during handling.

CLEANING AND STERILIZATION

Inspect each device before every use. The CUSP CALIPERS and ANNULAR SIZER should not be used if there are visible signs of cracking, crazing, deformation or other damage that may affect function.

Cleaning

The CUSP CALIPERS and SIZERS are supplied NON-STERILE and must be cleaned thoroughly and sterilized prior to initial use and each re-use. Discard all packaging materials before cleaning and sterilizing the CUSP CALIPERS and SIZERS. It is recommended that the tools be reprocessed as soon as reasonably possible following each use. The following manual cleaning and disinfection method is recommended for the CUSP CALIPERS and SIZERS.

1. Rinse devices with a solution of warm tap water and a hospital approved surfactant for approximately two (2) minutes.
2. Scrub the devices with a soft-bristled cleaning brush to remove all grossly visible debris.
3. Place the devices in an ultrasonic bath containing enzymatic detergent prepared according to manufacturer's directions and ultrasonically clean for fifteen (15) minutes.
4. If necessary, scrub the accessories using a cleaning brush to remove any remaining visible debris.
5. Rinse devices with warm tap water for approximately two (2) minutes to remove detergent.
6. Dry and place CUSP CALIPERS in the sterilization tray, close lid and wrap for sterilization.

The sterilization of the CUSP CALIPERS and SIZERS should be performed according to one of the following recommended steam sterilization methods:

Sterilization Parameters

Method (Steam)	Sterilizer Type	Minimum Temperature	Maximum Temperature	Minimum Time
Method 1 Preferred Method	Gravity Displacement	250°F (121°C)	254°F (123°C)	30 minutes
Method 2 Alternative method when immediate sterilization is required	Gravity Displacement	270°F (132°C)	275°F (135°C)	10 minutes
	Pre-Vacuum	270°F (132°C)	275°F (135°C)	4 minutes

Increasing the recommended sterilization temperature or time may result in reduced device life cycles.

Do not use aromatic solvent-based cleaners as disinfectant agents as these solvents may damage the components of the CUSP CALIPER or SIZER sets.

ORDERING INFORMATION

Model Number	Item
ART-2	Aortic Repair Tool Kit
ASR- 19/21	Annular Sizer size 19/21
ASR-23/25	Annular Sizer size 23/25
ASR-27/29	Annular Sizer size 27/29
ASR-31/33	Annular Sizer size 31/33
ARC-10	Cusp Caliper size 10
ARC-9	Cusp Caliper size 9
ARC-8	Cusp Caliper size 8
ARC-7	Cusp Caliper size 7
ARC-6	Cusp Caliper size 6
PST-ST5	Sterilization Tray

REFERENCES

Kunihara T, Aicher D, Rodionychewa S, Groesdonk H, Langer F, Sata F, Schäfers HJ. Preoperative aortic root geometry and postoperative cusp configuration primarily determine long-term outcome after valve-preserving aortic root repair. J Thorac Cardiovasc Surg 2012;143:1389-95

Bierbach B, Aicher D, Issa O, Bomberg H, Graber S, Glombitza P, Schäfers HJ. Aortic root and cusp configuration determine aortic valve function. Eur J Cardiothorac Surg 2010; 38: 400-406

Schäfers HJ, Schmied W, Psych D, Marom G, Aicher D. Cusp height in aortic valves. J Thorac Cardiovasc Surg. 2012;1-6

Svensson, Sizing for Modified David's Reimplantation Procedure. Ann Thorac Surg 2003;76:1751-3

Aicher D, Schneider U, Schmied W, Psych D, Kuniyara T, Tochii M, Schaefer HJ. Early results with annular support in reconstruction of the bicuspid aortic valve. J Thorac Cardiovasc Surg 2013;145 S30-4

Schäfers HJ, Bierbach B, Aicher D. A New approach to the assessment of aortic cusp geometry. J Thorac Cardiovasc Surg 2006;132:436-438

Lansac E, Di Centa I, Sleilaty G, Crozat E, et al. An aortic ring: From physiologic reconstruction of the root to a standardized approach for aortic valve repair J Thorac Cardiovasc Surg 2010;140:28-35

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO UNITED STATES CUSTOMERS ONLY:

ALTHOUGH THE AORTIC VALVE REPAIR TOOL KIT AND ACCESSORIES, HEREAFTER REFERRED TO AS "PRODUCT", HAVE BEEN DESIGNED, MANUFACTURED AND TESTED UNDER CAREFULLY CONTROLLED CONDITIONS, DUE TO CONDITIONS AND FACTORS THAT ARE OUTSIDE OF GBI'S SUPERVISION AND CONTROL AFTER SALE OF THE PRODUCT, GBI HEREBY DISCLAIMS ALL WARRANTIES EXPRESSED OR IMPLIED WITH RESPECT TO THE PRODUCT, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. GBI SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY DIRECT, INDIRECT, INCIDENTAL, CONSEQUENTIAL, SPECIAL OR OTHER TYPE OF DAMAGES INCLUDING, WITHOUT LIMITATION, ANY EXPENSES CAUSED BY OR RELATED TO ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER BASED UPON A WARRANTY, CONTRACT, TORT OR OTHER LEGAL THEORY. GBI NEITHER ASSUMES, NOR AUTHORIZES ANY OTHER PARTY TO ASSUME FOR IT, ANY WARRANTY, REPRESENTATION OR OTHER LIABILITY IN CONNECTION WITH THE PRODUCT OR ITS USE.

If any term or portion of this Disclaimer of Warranty and Liability provision is found to be unenforceable, illegal or contrary to applicable law by a court of competent jurisdiction, the remaining portions of this provision shall not be effected and shall remain valid, and any such unenforceability or illegality in jurisdiction shall not invalidate or render unenforceable this Disclaimer of Warranty and Liability provision in any other jurisdiction.

ENGLISH

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

DISCLAIMER OF WARRANTY ALTHOUGH THE AORTIC VALVE REPAIR TOOL KIT AND ACCESSORIES, HEREAFTER REFERRED TO AS "PRODUCT" HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE

PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. GBI, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. GBI SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

ENGLISH

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

Manufacturer:

Genesee BioMedical, Inc.
700 W. Mississippi Ave. Unit D5
Denver, CO 80223-4509, USA

Toll-Free: (800) 786 4890

Phone: (303) 777 3000

Fax: (303) 777-8866

www.geneseebiomedical.com

EU Authorised Representative:

Quality First International OÜ
Laki 30
12915 Tallinn
Estonia

UK Responsible Person:

Quality First International
1 Cook's Road
Stratford
London
E15 2PW

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LITERATUR

Kunihara T, Aicher D, Rodionychewa S, Groesdonk H, Langer F, Sata F, Schäfers HJ. Preoperative aortic root geometry and postoperative cusp configuration primarily determine long-term outcome after valve-preserving aortic root repair. J Thorac Cardiovasc Surg 2012; 143:1389-95

Bierbach B, Aicher D, Issa O, Bomberg H, Graber S, Glombitza P, Schäfers HJ. Aortic root and cusp configuration determine aortic valve function. Eur J Cardiothorac Surg 2010; 38:400-406

Schäfers HJ, Schmied W, Psych D, Marom G, Aicher D. Cusp height in aortic valves. J Thorac Cardiovasc Surg. 2012; 1-6

Svensson, Sizing for Modified David's Reimplantation Procedure. Ann Thorac Surg 2003; 76:1751-3

Aicher D, Schneider U, Schmied W, Psych D, Kunihara T, Tochii M, Schaefer HJ. Early results with annular support in reconstruction of the bicuspid aortic valve. J Thorac Cardiovasc Surg 2013; 145 S30-4

Schäfers HJ, Bierbach B, Aicher D. A New approach to the assessment of aortic cusp geometry. J Thorac Cardiovasc Surg 2006; 132:436-438

Lansac E, Di Centa I, Sleilaty G, Crozat E, et al. An aortic ring: From physiologic reconstruction of the root to a standardized approach for aortic valve repair. J Thorac Cardiovasc Surg 2010; 140:28-35

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