


# ROBOT<sup>x</sup> Robotic Tool Kit



## Tools for Robotic-Assisted Cardiac Surgery



**Genesee**  
BIOMEDICAL, INC.

## Suture Hooks

- Specialized tools for suture manipulation
- Available with left or right hook



## Specialized Suction Tools

- Improved access for chest wall and hard to reach areas
- Available in large diameter and ball tip



## Knot Pushers

- Stainless steel
- Open or Closed Tip




## Stainless Steel Valve Sizers

- Stainless-steel
- Sizes 25, 27, 29, 31, 33



## Flexible Silicone Sizers



## Ordering Information

### Robotic-Assisted Tools

Catalog Number	Description
RSH-L	Suture Hook Left
RSH-R	Suture Hook Right
RTK-BTS	Ball Tip Suction Tube
RTK-LDS	Large Diameter Suction Tube
RTK-KPO	Knot Pusher Open
RTK-RKP	Knot Pusher Closed
RTK-VSS	Stainless Steel Valve Sizer (Set of 5)
FRBS	Flexible Silicone Sizers (set of 9)
RPS-1850	Thoracic Port with Obturator (18mm ID, 50mm Length)
RPS-1870	Thoracic Port with Obturator (18mm ID, 70mm Length)
RPS-1890	Thoracic Port with Obturator (18mm ID, 90mm Length)
RPS-18110	Thoracic Port with Obturator (18mm ID, 110mm Length)

## Thoracic Ports



Robotic Tool Kit: Suture Hooks, Suction Tool, Robotic Ports and valve sizers  
**Indications:** The SUTURE HOOKS, KNOT PUSHERS, PORTS, OBTURATORS, SIZERS and SUCTION TOOLS are intended to aid in robotic-assisted cardiac surgery.  
 Use (IFU) for a full list of warnings, precautions and adverse events.  
 CAUTION: Federal law (USA) restricts this device to sale by, or on the order of, a physician.



Caution: Federal law (USA) restricts the sale, distribution, or use of this device to, by or on the order of a physician.  
 © 2022 Genesee BioMedical, Inc.  
 622001R00  
 All Rights Reserved.

[www.geneseebiomedical.com](http://www.geneseebiomedical.com)

For Customer Service:  
[cservice@geneseebiomedical.com](mailto:cservice@geneseebiomedical.com)

Phone: 303-777-3000  
 Toll free: 800-786-4890  
 Fax: 303-777-8866  
 Genesee BioMedical, Inc.  
 700 W Mississippi Ave, Unit D-5 Denver, CO 80223