

TruForm™ Sievers Annuloplasty Ring



Genesee
BIOMEDICAL, INC.

Shaped by Surgeons.
Built for Better Outcomes.

TruForm™ Sievers Annuloplasty Ring

TruePhysiologicalRepair.Semi-Rigid,DesignedforConsistencyAcrossthe Cardiac Cycle.

Designed For:



- ✓ Surgical reconstruction or remodeling of the mitral valve
- ✓ Procedures requiring durable, semi-rigid annular support
- ✓ Surgeons who prefer a ring congruent with surgical landmarks
- ✓ Anatomies involving the trigones, LVOT, and aortic annulus

Clinical Considerations

- ✓ Supports mechanical stress reduction through annular stabilization
- ✓ Aids in preserving physiologic motion with a semi-rigid design
- ✓ Intended for cases requiring precise anatomical fit across adjacent structures

Key Features:

- ✓ Designed to maintain optimum conditions throughout the cardiac cycle
- ✓ Reduces mechanical stress on the repaired valve
- ✓ Anatomically contoured to fit surgically relevant structures
- ✓ Reflects the shape of the mitral annulus and adjacent anatomy
- ✓ Braided polyester cover allows smooth suture passage
- ✓ Radiopaque markers assist with intraoperative orientation



Semi-Rigid Design. Anatomically Fit.

TruForm™ offers a semi-rigid ring that mirrors the natural shape of the mitral annulus while aligning with surgical landmarks for effective implantation and physiological support.

TruForm™ – Ordering Reference

Ordering Tip: Order using model number TFR-##

Example: TFR-30 for a full ring, size 30 mm

Catalog #	Product Description
TRH-24	TruForm™ Annuloplasty Ring size 24
TRH-26	TruForm™ Annuloplasty Ring size 26
TRH-28	TruForm™ Annuloplasty Ring size 28
TRH-30	TruForm™ Annuloplasty Ring size 30
TRH-32	TruForm™ Annuloplasty Ring size 32
TRH-34	TruForm™ Annuloplasty Ring size 34
TRH-36	TruForm™ Annuloplasty Ring size 36
TRH-38	TruForm™ Annuloplasty Ring size 38
TRH-40	TruForm™ Annuloplasty Ring size 40
TRH-42	TruForm™ Annuloplasty Ring size 42

Stainless-Steel
Handle



Stainless-Steel
Sizer



TruForm™ – Ordering Reference

Ordering Tip: Order using the model number followed by the size.

Example: TRH-32 for a TruForm™ Ring, size 32 mm

Catalog #	Accessory Description
FRBS	Flexible Silicone Sizers, set of 9 (Sizes 24–40)
GRH-U	Reusable Handle
GRH-UXL	Extra-Long Reusable Handle (35 mm)
GSS-26/28	Stainless Steel Sizer, double-ended 26/28
GSS-30/32	Stainless Steel Sizer, double-ended 30/32
GSS-34/36	Stainless Steel Sizer, double-ended 34/36
GSS-38/40	Stainless Steel Sizer, double-ended 38/40
GSS-4	Stainless Steel Sizer, set of 4

Indications: The TruForm™ Annuloplasty Rings are indicated for the surgical reconstruction or remodeling of diseased or damaged mitral valves. The ring provides support for and restricts expansion of the mitral annulus.
 Warnings: Only surgeons who have received adequate training to determine whether incompetent, stenotic, or diseased heart valves are capable of being repaired or replaced should use this device. Refer to the Instructions for Use (IFU) for a full list of warnings, precautions and adverse events.
 CAUTION: Federal law (USA) restricts this device to sale by, or on the order of, a physician.
 Image Credits: Images courtesy of Dr. Kevin Miller, Dr. Steven Bolling, and Dr. Jose Navia.





TruForm™ Sievers Annuloplasty Ring | Semi-Rigid Mitral Valve Repair with An

Genesee Biomedical



Genesee
BIOMEDICAL, INC.

TRUFORM™ SIEVERS

MITRAL ANNULOPLASTY



SUBSCRIBE



Watch on YouTube



Annuloplasty Accessories Quick Reference

Simplify. Standardize. Support.

SKU Reference Table

Catalog #	Accessory Description	Compatible Products
FRBS	Set of 9 Flexible Silicone Sizers (Sizes 24–40)	FlexForm™, WellsForm™, TruForm™
GSC	Annuloplasty Scale	TransForm™, FlexForm™, WellsForm™, TruForm™
GSC-2	Set of 2 GSC Annuloplasty Scales and 2 handles (1 GRH and 1 GRH-U)	TransForm™, FlexForm™, WellsForm™, TruForm™
GSS-26/28	Stainless Steel Sizer, double-ended 26/28	FlexForm™, WellsForm™, TruForm™
GSS-30/32	Stainless Steel Sizer, double-ended 30/32	FlexForm™, WellsForm™, TruForm™
GSS-34/36	Stainless Steel Sizer, double-ended 34/36	FlexForm™, WellsForm™, TruForm™
GSS-38/40	Stainless Steel Sizer, double-ended 38/40	FlexForm™, WellsForm™, TruForm™
GSS-4	4 Stainless-Steel Sizers (26/28, 30/32, 34/36, and 38/40) with 2 GRH-U handles	FlexForm™, WellsForm™, TruForm™
GRH	Tricuspid Stainless Steel Sizers, double-ended	FlexForm™, WellsForm™
GRH-U	Reusable Handle	TransForm™, FlexForm™, WellsForm™, TruForm™
GRH-UXL	Extra-Long Handle (35 cm)	TransForm™, FlexForm™, WellsForm™, TruForm™
TFRS	Set of Reusable Silicone Sizers	TransForm™
TFSS	Set of Stainless-Steel Sizers (includes 1 GRH and 1 GRH-U)	TransForm™
WFBS	Flexible Silicone Sizers, set of 6, Sizes 26-36	FlexForm™, WellsForm™

Genesee BioMedical offers reusable sizing tools and handles to support precise valve repair across our annuloplasty portfolio. While most accessories are shared among FlexForm™, WellsForm™, and TruForm™, **TransForm™ utilizes its own sizers** to ensure accurate mitral remodeling.*

Refer to individual product sections or Instructions for Use (IFU) for full cleaning and compatibility details.

All Genesee annuloplasty accessories featured in this section are fully **reusable** and designed for repeat sterilization

* Accessories must be resterilized before each use. Refer to the product's Instructions for Use (IFU) for full warnings, precautions, and reprocessing instructions.
CAUTION: Federal law (USA) restricts these devices to sale by, or on the order of, a physician.

Contact

Shaped by Surgeons.
Built for Better Outcomes.



Genesee BioMedical, Inc.
700 W. Mississippi Avenue, Unit D5
Denver, CO 80223-4509



Phone:
303-777-3000



Website:
www.geneseebiomedical.com



General Inquiries:
info@geneseebiomedical.com

“Innovation begins with listening. Every great product we build starts with a surgeon’s challenge.”
— Matt Dambeck, CEO Genesee BioMedical, Inc

